



THE AMERICAS

Owen Mumford Inc.
1755 West Oak Commons Court
Marietta, Georgia 30062, USA
T : (800) 421-6936 F : (770) 977-2866
E : info@owenmumfordinc.com

**HEAD OFFICE
INTERNATIONAL SALES & MARKETING**

Owen Mumford Ltd.
Brook Hill, Woodstock
Oxfordshire OX20 1TU, United Kingdom
T : 01993 812021 F : 01993 813466
E : info@owenmumford.co.uk

UK / IRELAND SALES & MARKETING

T : 01608 645555 F : 01608 646546
E : info@owenmumford.co.uk

FRANCE

T : 02 32 51 88 70 F : 02 32 21 96 95
E : contact@owenmumford-fr.com

GERMANY

T : 06026/97750 F : 06026/977529
E : info@owen-mumford.de

CHINA

T : 021 6032 3528 F : 021 6032 3782
E : info@owenmumford.co.uk

www.owenmumfordinc.com


OUR ACCREDITATIONS INCLUDE:

- ISO13485 Certification
- Medical Device Directive 93/42/EEC
- CMDCAS 13485 (Canadian Health Authority)
- Environmental Management Systems ISO14001




On-the-Go® Use Has Never Been Easier

Unifine® Pentips® Plus are easy to use, convenient and safer† than most pen needles when traveling or injecting On-the-Go®.



Remove protective seal, screw on NEW pen needle and discard needle shield.




After injection, push USED pen needle into removal chamber until it clicks, then unscrew and dispose of safely.



Scan for more information and instructional video.

Selecting an Injection Site

Where you inject on the body will affect the rate of drug absorption. Rotating the site you use will make injections easier, safer and more comfortable. Repeated use of the same site can cause hard areas under the skin preventing the insulin from being properly absorbed. Use the diagrams provided to help determine the best injection sites for you, and to track where you are injecting.



Users should seek medical advice regarding injection sites, technique or choice of pen needle.

Range and Compatibility



Unifine® Pentips® Plus

100ct 08470-3840-01 100ct 08470-3850-01 100ct 08470-3890-01 100ct 08470-3830-01 100ct 08470-3829-01

Fits ALL Leading Brands of Injection Pens.



All trademarks are the property of their respective owners.



A Safer†, More Convenient Pen Needle

Unifine® Pentips® Plus



The Only Pen Needle with...

Built-In Remover



For further information visit: www.unifinepentipsplus.com

UPPL17/OM/0114/1/US

Unifine[®] Pentips[®] Plus

The Next Generation in Pen Needles

Unifine[®] Pentips[®] Plus deliver the ultimate in convenience combining the superior performance of Unifine[®] Pentips[®] with a built-in remover for safer On-the-Go[®] use.

Ergonomic Design

Unifine[®] Pentips[®] Plus makes pen needle attachment and removal easy, with a larger container for improved grip and control.

Easy Removal for On-the-Go[®]

Click used needles into locking chamber for hassle-free removal, designed to reduce the risk of accidental needlestick.



Disposal Chamber

When sharps containers are not convenient, Unifine[®] Pentips[®] Plus offers a built-in locking chamber designed to hold used needles until they can be disposed of safely. Exclusive *Safety Click Technology*[®] clicks to indicate the used pen needle is locked in the chamber.

[†]Users agreed that Unifine Pentips Plus helped reduce risk of needlestick injury versus their prior pen needle. Data on file.

Removing used pen needles after every injection is **important**. The unique design of Unifine[®] Pentips[®] Plus makes removing needles easy and convenient.

Reuse of pen needles can cause complications:

- Painful injections, bruising and Lipohypertrophy
- Insulin crystallization and blockage
- Air bubbles and drug leaking in and out of the cartridge, affecting dose accuracy



Drug Delivery

Drug Delivery



Healthcare Professional Preferred

Diabetes nurses in our research thought Unifine[®] Pentips[®] Plus provided safer removal and storage of used pen needles than standard pen needles.

Data on file

Exceeding Expectations for Quality and Performance

For over 60 years, Owen Mumford has established a global reputation for innovation and excellence in medical device manufacturing. Produced under strict quality standards, Unifine Pentips Plus undergo 100% inspection to ensure our pen needles offer the highest quality and performance.

The Unifine[®] Pentips[®] Advantage

- Triple sharpened for one of the **lowest penetration forces**
- Anti-coring technology for **improved drug delivery**
- Advanced silicone for a **lower glide force**
- Thin wall technology for **improved drug flow** and **low delivery force**





"I find Unifine[®] Pentips[®] Plus easy to use and much safer. You don't have to get near the needle when you change it and insert it back into the holder. It's safe without having to touch it at all."


Male, diagnosed 14 years ago, injects 3 times a day

Proven Performance

New Unifine[®] Pentips[®] Plus represent the next generation in pen needle technology, combining the superior performance of Unifine[®] Pentips[®], with the world's first built-in remover for safer, personal use. Whether On-the-Go[®] or using at home, Unifine[®] Pentips[®] Plus deliver exceptional quality and convenience.



Unifine[®] Pentips[®] Plus With Anti-Coring



Other Pen Needle Without Anti-Coring

Tried and Tested

At Owen Mumford we don't just test our products in the lab, we talk to patients and healthcare professionals too. In a recent study **participants preferred Unifine[®] Pentips[®] Plus to their current needle[†]**, describing Unifine[®] Pentips[®] Plus as a safer, easy-to-use pen needle that reduced the risk of needlestick injury.

[†]Data on file.

