# Flu Season is Almost Here: Are You Prepared?

Prevent unnecessary delays and costs by stocking puritan UniTranz-RT® transport systems for flu season.

All Puritan brand products are made in the USA, giving you more reliable availability!



# PURITAN UNITRANZ-RT® & FLOCKED SWABS:

Safe, effective, and easy to use, Puritan's UniTranz-RT° transport systems are versatile in culturing viruses, chlamydia, mycoplasma, and ureaplasma samples.

- Medium remains stable when stored at room temperature
- Maintains organism viability at room temperature
- Self-centering, leak-proof cap
- Fully compatible with automation systems, EIA, PCR, DFA and molecular assays
- · Easy swab removal & disposal
- Various swab-handle breakpoints
- Available in 1ml & 3ml volumes



Find out about pricing and availability by contacting our sales representatives today:

+1 207-876-3311

PURITAN UNITRANZ-RT° TRANSPORT SYSTEM (3ML FILL & ONE ELONGATED FLOCKED SWAB) >



Puritan's UniTranz-RT°, a universal transport medium system, was formerly known as "UTM-RT." To better reflect the quality of our distinct products, Puritan has renamed this product line as UniTranz-RT°. Puritan's UniTranz-RT° product line has not changed and new validation is not necessary.






UT-306

# PurFlock Ultra® Swabs Mean Better Absorption and Elution

Puritan's UniTranz-RT° transport systems for flu testing utilize proprietary PurFlock Ultra° swabs. Flocked swabs are recommended by the American Society for Microbiology and are suitable for RIDT and viral cell culture.



# PURITAN UNITRANZ-RT° TRANSPORT SYSTEM WITH UNIVERSAL TRANSPORT MEDIUM

Product	Vial Length	Medium Volume
UT - 100	12 x 80 mm	1ml
UT - 300	16 x 100 mm	3ml



Our knowledgeable sales representatives are standing by to answer any questions you might have about pricing and availability:

+1 207-876-3311

## ORDERING INFORMATION

### PURITAN UNITRANZ-RT° TRANSPORT SYSTEM WITH PATENTED PURFLOCK ULTRA° SWABS AND UNIVERSAL TRANSPORT MEDIUM

Product	Vial Length	Medium Volume	Swab Type (Sterile)	Breakpoint
UT - 106	12 x 80 mm	1ml	1 large PurFlock Ultra®	80 mm
UT - 116	12 x 80 mm	1ml	1 mini-tip PurFlock Ultra®	80 mm
UT - 117	12 x 80 mm	1ml	1 ultrafine PurFlock Ultra®	100 mm
UT - 306	16 x 100 mm	3ml	1 elongated PurFlock Ultra®	100 mm
UT - 316	16 x 100 mm	3ml	1 mini-tip PurFlock Ultra®	100 mm
UT - 317	16 x 100 mm	3ml	1 ultrafine PurFlock Ultra®	100 mm
UT - 367	16 x 100 mm	3ml	1 ultrafine PurFlock Ultra®, 1 elongated PurFlock Ultra®	100 mm

# PURITAN UNITRANZ-RT° TRANSPORT SYSTEM WITH SPUN POLYESTER SWABS AND UNIVERSAL TRANSPORT MEDIUM

Product	Vial Length	Medium Volume	Swab Type (Sterile)	Breakpoint
UT - 302	16 x 100 mm	3ml	1 regular polyester, 1 mini-tip polyester	100 mm
UT - 361	16 x 100 mm	3ml	1 regular polyester	100 mm
UT - 362	16 x 100 mm	3ml	2 regular polyester	100 mm
UT - 366	16 x 100 mm	3ml	1 regular polyester, 1 mini-tip polyester	100 mm



